

SpotOptics

The software people for optics

LENTINO

AUTOMATED WAVEFRONT SENSOR

- **Accurate Metrology of standard and aspherical lenses**
- **$\phi=0.3$ to $\phi=20$ mm**
- **F/1 to F/15**
- **Accurate motor for z-movement**
- **Accurate XY and tilt stages for easy centering of lenses**



LENTINO

TECHNICAL SPECIFICATIONS

<i>Focal ratios covered</i>	f/1 to f/15 standard
<i>Diameter of lens that can be tested in parallel light</i>	20mm
<i>Diameter of lens that can be tested with pinhole at focus</i>	15mm
<i>No of spots (with standard camera)</i>	23x23
<i>Maximum no. of spots (large format camera)</i>	80x80
<i>Diameter and focal length of standard lenslets</i>	(0.2,11),(0.2,22),(0.3,41) mm
<i>Measure aspherical elements</i>	~15% variation in longitudinal spherical aberration
<i>RMS repeatability of computation of Zernike polynomials</i>	1-2nm rms ($\lambda/600$ - $\lambda/300$) at 633nm
<i>RMS repeatability of wavefront measurements</i>	$< \lambda/200$
<i>Accuracy</i>	$\lambda/10$ - $\lambda/100$ – depending on the accuracy of the calibration elements
<i>Dynamic range of measurements sub-pupil (tilt subtracted)</i>	$\pm 50 \lambda$
<i>Wavelength range</i>	UV (0.193-1.1 μ), Vis (0.193-1.1 μ), 0.95-1.7 μ
<i>Light sources</i>	LEDs, LDs and Halogen lamps available (at one wavelength or in white light)
<i>Software (control and analysis)</i>	Sensoft
<i>Camera - connections</i>	Cameras with Gigabit and USB2 and USB3 connection available. 10-bits to 16-bits
<i>Acquisition speed</i>	15-500Hz (camera dependent)
<i>Processing speed</i>	5-150 Hz (camera and PC dependent)
<i>Power requirement for stepper motor</i>	24V/2A DC



SENSOFT: THE SOFTWARE

Sensoft: The modular software package:


- Fully controls the hardware of Lentino
- Performs the Shack-Hartmann (SH) analysis
- Computes Zernike coefficients, diagnostics (alignment and correct focal plane), wavefront, MTF, spot diagram
- Has a Loop mode for on-line adjustment of optical system

Lentino in your production line:

- Lentino – with its own PC - can easily be adapted to the production line
- It can work in a closed-loop with the PC of the manufacturing machine
- A software module defines the communication protocol and transfers the results between the PCs




ON-LINE ALIGNMENT OF COMPLEX OPTICAL SYSTEM IN A FAST LOOP



Alignment started
Coma decreases

- The alignment of complex optical systems becomes easy by monitoring coma and astigmatism in a continuous loop
- The individual (x, y) components of coma and astigmatism, as well as the total coefficients are displayed
- The optimization can be done for one component at a time, as the software can display one component of interest
- Optimal alignment is reached when the coma and astigmatism components converge towards a given tolerance


MTF MEASUREMENTS



MTF after subtracting the contributions of tilt and defocus present in the data

LENTINO

LAYOUT



PHYSICAL

Dimensions: 20 (L) x 29 (W) x 45 (H) cm
 Weight: ~5 Kg
 Cameras: USB3, Gigabit Ethernet
 Motor power supply: 24V, 2A

KEY FEATURES

- Measurement technique**
 Shack-Hartmann wavefront sensor
- Test in parallel light or at the lens focus**
 Parallel light (with a collimator)
 At the focus of the lens (with a pinhole)
 Light sources with different wavelength available
- Calibration**
 High-quality parallel light source
 Pinhole calibration unit
- Easy access to optical elements**
 Easy access to optical group for on-line alignment

SOFTWARE

Easy alignment of lens group via software: the software gives graphical indication of the misalignment of the optical system, using coma and astigmatism.

May 2014

www.spot-optics.com