# Falcon III EMCCD - NEW

Digital Monochrome Scientific Frame Transfer EMCCD  $1024 \times 1024 \cdot 10 \mu m \times 10 \mu m$  pixels • Cooled to  $-100^{\circ}$ C • 1MP Scientific •





### **Key Features and Benefits**

NEXT GENERATION photon counting sensitivity

- Lower read noise of <0.01e-Best sensitivity of any camera technology
- Faster readout in full resolution
   x 3 times faster than previous generations
- Higher EM gain of x 5000
   See single photon events
- Up to 95% QE from back-illuminated sensor Optimum Photon collection
- Strong UV and NIR reponse and ultrawide bandwidth From 200nm through to 1100nm
- Deep cooled to -100°C
   For minimal background events

### EMCCD - GEN III A NEW GENERATION

## The Photon Harvester!

Resolution	1024 x 1024
Pixel Size	10μm x 10μm
Readout Noise	<0.01e-
Frame Rate	34fps
Cameralink	16bit



#### **Specification for Falcon III EMCCD**

Sensor Type	1" Back Thinned Frame Transfer EMCCD
Active Pixel	1024 x 1024
Pixel Size	10µm x 10µm
Active Area	10.2mm x 10.2mm
Full Well Capacity	35,000 electrons
Shift Register Well Depth	200,000 electrons
Non-Linearity	<1%
Readout Noise	<0.01 electrons with EM gain ON,
	<60 electrons with EM gain OFF
Full Resolution Frame Rate	34fps
Dark Current (@ -100°C)	0.0002 e/p/s
Digital Output Format	16 bit CameraLink (base configuration)
Peak Quantum Efficiency	95%
Spectral Response	200 - 1100nm
Cooling	-100°C with +10°C coolant
Binning	1x1 up to 32x32
Lens Mount	C-Mount
Synchronisation	Trigger IN and OUT - TTL compatible
Power Supply	12V DC ±10%
Total power consumption	<100W
Operating case temperature	-20°C to +55°C
Storage Temperature	-30°C to +85°C
Dimensions	129mm x 112mm x 94mm
Weight (no lens)	<1.5Kg

Raptor Photonics Limited reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.

#### **Ordering Information**

#### Camera

Falcon EM351 digital B/W camera RPL-FA351V-BV-CL Falcon Power Supply Cable RPL-KY-CBL

#### **Optional Accessories**

EPIX(R) EB1 base CL card RPL-EPIX-EB1

EPIX(R) base Notebook CL card RPL-EPIX-ECB1-34

EPIX(R) base Notebook CL card RPL-EPIX-ECB1-54

EPIX(R) XCAP STD software RPL-XCAP-STD

CameraLink Cable, 2m<sup>2</sup> RPL-CL-CBL-2M

Optical Visible lenses<sup>3</sup> RPL-xx-xxxx

Note 1: Extended operating temperature range on request.


Note 2: Longer CL cable available up to 25M

Note 3: Please consult us to check our range of lenses

Demo is available on request. Pricing AOR subject to volumes.

Detailed technical drawings can be downloaded at www.raptorphotonics.com

#### **Quantum Efficiency**



#### **Applications**

- Adaptive Optics and Astronomy
- Calcium signaling
- Fluorescence imaging / spectroscopy
- Flow cytometry
- FRET / FRAP / TIRF
- · Genome sequencing
- · High content screening
- · High resolution fluorescence imaging
- Hyperspectral imaging
- · Live cell imaging
- Photon counting
- Single molecule detection
- Solar cell inspection
- X-ray & High energy



Willowbank Business Park Larne, Co Antrim BT40 2SF, Northern Ireland ROW Sales T: +44(0)2828 270 141 E: sales@raptorphotonics.com www.raptorphotonics.com USA Sales T: +1 (770) 364-7240 E: request@phxatl.com www.phxatl.com



Document #: INFA351V-BV-CL 0117R1