Fiber Optic Raman Probe






Ball-lens Raman fiber-optic probe is a new addition to art photonics' FlexiSpec® product line.

The technology of direct deposition of the dielectric filters at the fiber end faces enables miniaturization of the fiber-optic probe head without compromising its performance. It makes possible to develop a cost-effective Raman probe for endoscopic and other applications.

The fiber optic Raman probe is produced for multiwave excitation in the range 690-785nm, e.g. @785nm – "Finger print" spectral range with fluorescence reduction, and @690nm – "High wavenumber" spectral range for conventional Raman spectrometers.

Applications:

- Biomedical
- Biotechnology
- Material Analysis
- Food and Drink
- Chemical
- Pharmaceutical
- Petrochemical
- R & D



1x 9 Fiber bundle Dielectric filters for excitation and detection channels directly deposited on fiber ends.



1x Excitation Fiber NIR200/ 220/AL/FC/PC



9x Detection Fiber NIR200/ 220/AL SMA905 (line bundle)



Specification	
Laser Wavelength:	690-785nm
Available Spectral Range:	2000-3800 cm ⁻¹ (690nm), 200-2000 cm ⁻¹ (785nm),
Fibers	1+ 9 (one excitation & 9 collection fibers with 200µm core dia.)
Long-pass filter transmission range	810-950nm (>90%)
Long-pass filter blocking range	690-800nm (OD>6)
Short-pass filter transmission range	690-790nm (>90%)
Short-pass filter blocking range:	810-950nm (OD>6)
Ball Lens	2mm diameter silica
Total Length	1.5m
Shaft (needle) Length	150mm
Shaft (needle) Diameter	3mm
Length of Legs	0.3m
Minimal Bending Radius	130mm
Input / Output Connectors	Excitation – FC/PC; Detection – SMA905
Shaft	Stainless Steel
Protective tubing	Silicon coated Stainless Steel Conduit
Temperature range	-40 to 150°C
Pressure limits	0-7 Bar (100 PSI)

